HC220^{SERIES}

HC220LE/HC221LE

with Ambient Tracking™ PLUS

English-Francais-Español

CE 0123

Patents Pending in Principal Countries of the World

Fisher & Paykel
HEALTHCARE

HC220^{SERIES}

HC220LE/HC221LE

SECTION

English · · · · · A

Francais · · · · · B


Español · · · · · C

HC220 SERIES

HUMIDIFIED CPAP SYSTEM

with *Ambient Tracking™ PLUS*

English

The HC220 Series (HC220LE / HC221LE) humidified CPAP system is an integrated CPAP and heated humidifier with *Ambient Tracking™ PLUS* technology for use in the home.

CAUTION: USA FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

CONTENTS

| | |
|---|---|
| Symbol Definitions | 2 |
| Warnings | 2 |
| Cleaning & Maintenance | 2 |
| Grounding Instructions | 3 |
| Package Contents | 3 |
| Important Parts of Your HC220 Series System | 3 |
| Getting Started | 4 |
| Making the Most of Your HC220 Series System | 4 |
| Using the HC220 Series Warm-up Mode | 4 |
| Controlling Humidity | 5 |
| <i>Ambient Tracking™ PLUS</i> | 5 |
| Compliance Monitoring (HC221LE only) | 5 |
| Frequently Asked Questions | 6 |
| Product Specifications | 6 |
| Troubleshooting | 6 |

Symbol Definitions

Warm-up
Mode

Type B

ATTENTION
Consult
accompanying
documents

Alternating
Current

IPX1
Drip Proof

CE 0123

CE Mark
93/42/EEC
Class IIb

Humidity
Control

Double
Insulated

WARNINGS

READ ALL INSTRUCTIONS BEFORE USE

DANGER - To reduce the risk of electrocution:

- Always disconnect the HC220 Series immediately after using.
- Do not use the HC220 Series while bathing.
- Do not store the HC220 Series where it can fall or be pulled into water.
- Do not place or drop the HC220 Series into water or other liquid.
- Do not reach for an HC220 Series that has fallen into the water. Unplug immediately.

WARNING - To reduce the risk of burns, electrocution, fire or injury:

- Do not touch the heater plate. The heater plate temperature may exceed 65° C(149° F)
- The HC220 Series should not be left unattended when plugged in.
- Close supervision is necessary if the HC220 Series is used by, on, or near children or invalids.
- Use the HC220 Series only for its intended use as described in this manual.
- Never operate the HC220 Series if it has a damaged power cord or plug, if it is not working properly, if it has been dropped or damaged, or any part of the HC220 Series has been dropped into water. Return the HC220 Series to an authorised service center for examination and repair.
- Keep the power cord away from heated surfaces.
- Never drop or insert any object into any opening or tube.
- Do not use outdoors or operate where aerosol (spray) products are being used.
- Connect the HC220 Series to a properly grounded outlet only. See *Grounding Instructions*.
- **Remove humidification chamber before filling, and fill humidification chamber with distilled water only.** See *Getting Started*.
- At low CPAP pressures, the flow through the mask may be insufficient to clear all exhaled gas from the tubing, and consequently some rebreathing may occur.
- Adjusting the pressure control can result in ineffective therapy. The pressure control must only be adjusted by the clinician.
- Equipment not suitable for use in the presence of flammable anaesthetic, mixed with air or oxygen or nitrous oxide.

- The HC220 Series should only be used with masks and connectors recommended by Fisher & Paykel Healthcare, or by your physician or respiratory therapist. A mask should not be used unless the HC220 Series is switched on, and operating properly. The vent hole(s) associated with the mask should never be blocked.

CAUTION

- Do not use your HC220 Series in close proximity to high frequency apparatus, shortwave or microwave equipment. Emissions from such devices may effect the HC220 Series' operation.
- Never block the air openings of the HC220 Series or place it on a soft surface such as a bed or couch, where the openings may be blocked. Keep the air openings free of lint, hair or the like.

OXYGEN USE

If supplemental oxygen is used with the HC220 Series, the oxygen flow must be turned off when the HC220 Series is not operating. This ensures that oxygen will not flow through the tube and accumulate in the HC220 Series enclosure.

- The HC220 Series should be located in a position where ventilation around the device is not restricted.
- It is recommended that supplemental oxygen be administered at the mask.
- Oxygen should not be used while smoking or in the presence of an open flame.
- The use of oxygen in therapy requires that special care be taken to reduce the risk of fire. Any materials which will burn in air and some that will not are easily ignited and burn rapidly in high concentrations of oxygen. Accordingly, for safety it is necessary that all sources of ignition be kept away from the product and preferably out of the room in which it is being used. "NO SMOKING" signs should be prominently displayed.
- A spontaneous and violent ignition may occur if oil, grease or greasy substances come in contact with oxygen under pressure. These substances must be kept away from oxygen regulators, cylinder valves, tubing and connections, and all other oxygen equipment.
- At a fixed flow rate of supplemental oxygen flow, the inhaled oxygen concentration will vary, depending on the pressure settings, patient breathing pattern, mask selection and leak rate.

DO NOT DISCARD

Cleaning & Maintenance

NOTE: The HC220 Series should be cleaned as required.

1. Switch off the HC220 Series and unplug from the power socket.
2. Wipe the exterior of the HC220 Series with a clean, damp (not wet) cloth and mild dishwashing detergent. Note: do not use harsh abrasives or solvents, or these may damage the HC220 Series.
3. Daily - Rinse the water chamber and breathing tube with warm, soapy water.
4. Weekly - Soak the inside of the chamber for 10 minutes in a solution of one part white vinegar to ten parts water.
5. Replace the air filter when it becomes significantly discolored.

WARNING

- Refer all other maintenance to qualified service personnel.
- Always switch off the HC220 Series and unplug the power cord from the power socket before cleaning.
- Never rinse or place the unit in water or allow water to get into the interior of the unit.
- Ensure that the unit is dry before plugging into the power socket.

Grounding Instructions (for CLASS I models only)

NOTE: Models labeled as Class II (Double Insulated) are marked with the  symbol and do not require grounding. The following applies to Class I (three wire) models only.

The HC220 Series should be grounded. In the event of an electrical short circuit, grounding reduces the risk of electric shock by providing an escape wire for the electric current. The HC220 Series is equipped with a cord having a grounding wire with a grounding plug. The plug must be plugged into an outlet that is properly installed and grounded.

DANGER - Improper use of the grounding plug can result in a risk of electric shock.

Check with a qualified electrician or serviceman if the grounding instructions are not completely understood, or if in doubt as to whether the HC220 Series is properly grounded.

Extension Cords:

If it is necessary to use an extension cord, use only a 3-wire extension cord that has a three-blade grounding plug, and a 3-slot receptacle that fits the HC220 Series plug. Replace or repair a damaged cord.

In Canada/USA only:

This product is for use on a nominal 120V circuit, and has a grounding plug that looks like the plug illustrated (Figure 1). A temporary adaptor, may be used to connect this plug to a 2-pole receptacle if a properly grounded outlet is not available. The temporary adaptor should be used only until a properly grounded outlet (Figure 1) can be installed by a qualified electrician. The green coloured rigid ear, lug, etc extending from the adaptor, must be connected to a permanent ground such as a properly grounded outlet box cover. Whenever the adaptor is used, it must be held in place by the screw.




Figure 1: Grounding methods

If repair or replacement of the cord or plug is necessary, do not connect the grounding wire to either flat blade terminal. The wire with insulation having an outer surface that is green with or without yellow stripes is the grounding wire.


Package Contents

Accessories

Important Parts of Your HC220 Series System

Getting Started

1. Remove your HC220 Series from its packaging.
2. Place the HC220 Series on a low shelf or on the floor beside your bed. Warning: Ensure that the HC220 Series is placed at a level lower than head height while sleeping.
3. Remove one water chamber from the packaging.
4. Remove the blue caps and discard. (See Fig. 1)
5. Using the funnel provided, fill the water chamber up to the black line with distilled water. (See Fig. 2)
NOTE: Never fill the water chamber while it is attached to the HC220 Series. When moving your HC220 Series, ensure that the water chamber is empty. Machine failure due to water damage is not covered by warranty.
6. To attach the water chamber to the HC220 Series, press down the finger guard and slide the chamber on. (See Fig. 3)
7. When the chamber is fitted correctly, the finger guard will click into place. (See Fig. 4)
8. One end of the tube has a white plastic connector. Push this onto the outlet on top of the chamber. (See Fig. 5)
9. Connect the other end of the tube to the mask.
10. Plug the HC220 Series power cord into your household power socket.
11. Set the humidity control on the front of the HC220 Series to '2' (or as recommended by your homecare provider or sleep lab). If heated humidification is not required, set the control to the **off region**.
12. Switch on the HC220 Series by pressing the power switch on the right hand side. (See Fig. 6)



Fig 1


Fig 2


Fig 3


Fig 4


Fig 5


Fig 6

The HC220 Series is now ready for use.

Making The Most Of Your HC220 Series System

Part of the Fisher & Paykel Healthcare philosophy is to design innovative and effective features to ensure your CPAP treatment is always comfortable and simple. The following unique features are available with your CPAP unit

USING THE HC220 SERIES WARM UP MODE


Some users, when first using CPAP therapy, find the airflow uncomfortable. To help you fall asleep comfortably, the HC220 Series provides a gradual increase in pressure and humidity if required. We call this feature **warm-up**. Warm-up can be used with or without heated humidification. To use warm-up, follow these simple steps:

1. When the HC220 Series is switched on (See 'Getting Started'), it will **automatically deliver your set pressure** (the green warm-up light will be **off**).


To activate **warm-up mode** simply press the warm-up button on the front of the HC220 Series (the green warm-up light will turn **on**).

When **warm-up mode** is activated, it will allow you to fall asleep as the unit gradually increases the pressure over a period of 20 minutes until set pressure is reached. Pressure will start at a minimum of 3cmH2O (if pressure is less than 9cmH2O) or 1/3 of set pressure.

Note: To deactivate warm-up mode, simply switch off the HC220 Series (See 'Getting Started'), pause for a few seconds, then switch it back on.


2. Activating warm-up mode during the night: If you want to activate **warm-up mode** during the night, simply press the warm-up button on the front of the HC220 Series.

Controlling Humidity

CPAP therapy often dries out the nasal passages, sometimes causing discomfort for the user. Clinical studies have shown that heated humidification helps to reduce this discomfort and can help the user to benefit more from CPAP therapy.

You are able to adjust the amount of humidity that the HC220 Series generates by adjusting the humidity control on the front of the unit. Occasionally, some users prefer not to use heated humidification. Heated humidity can be disabled by setting the humidity control to the 'off' region.

Ambient Tracking™ PLUS

Ambient Tracking™ PLUS maintains humidity due to changes in room temperature and airflow without condensation collecting in the tube. The heater plate automatically adjusts for changes in room temperature and air leaks caused by mouth leak and/or mask leak to maximize humidity and minimize condensation.


Compliance Monitoring

Sometimes CPAP users find that when they wake in the morning, they have taken their CPAP mask off while they were sleeping. CPAP therapy is only effective while the user is wearing the mask.

Clinical studies have shown that the use of heated humidification increases the amount of time that CPAP is used overnight, which is why Fisher & Paykel Healthcare have designed the HC220 Series.

Compliance monitoring provides you and your medical practitioner with the information you need to assess your progress with CPAP therapy.

When the HC220 Series is first turned on, the *compliance display* will show the following pieces of information for 8 seconds each before returning to the default display - set pressure.

1. Machine on time (identified by the **hr** symbol)

This is the total number of hours that the HC220 Series has been **switched on** (*this includes hours with the mask on and hours with the mask off*).

Example:

634 hr

This means that the HC220 Series has been switched on for a total of 634 hours.

2. Average hours of use per night (identified by the **hr/** symbol) (HC221LE only)

This is the average number of hours that the HC220 Series has been used with the mask on each night.

Example:

7.4 hr/ 

This means that the HC220 Series has been used with the mask on for an average of 7.4 hours per night.

3. Check sum (identified by the symbol) (HC221LE only)

This is a number which can help your medical practitioner to check that the other values are correct.

Example:

 564

Your medical practitioner or homecare provider may ask you to give them this number for reporting purposes.


4. Set pressure (identified by the **cm H₂O** symbol)

This is the pressure which the HC220 Series is set to in cm of H₂O – for pressure measurement.

Example:

9.5 cm H₂O

This means that the HC220 Series has been set to a pressure of 9.5 cm H₂O. Note: The display will continue to show only set pressure. To see the other information, simply turn the HC220 Series off, pause for 3 seconds, then turn it back on.

(NB) If the key symbol [] is displayed alongside set pressure, it has been locked and cannot be adjusted. Please consult with your medical practitioner or Homecare provider should you wish to change your recommended therapeutic pressure setting.

Frequently Asked Questions

When I wake up in the morning, my nose and throat feel dry. What can I do?

Try increasing the setting on your humidity control. If this does not help, please contact your clinician for advice.

How do I prevent condensation in the tubing?

The humidity control enables adjustment of humidity, so that there should be fine misting in the six inches of the breathing tube closest to your face. If the humidity control is set too high for the conditions, condensation may occur in the breathing tube. Note: Using greater than 6' (183 cm) of tubing will increase condensation.

There are several ways of reducing condensation in the tube:

1. Ensure the HC220 Series is not positioned in a cool draft.
2. Insulate the tubing by positioning it under the bedcovers, by using the 900HC003 rainout reduction kit, wrapping it with flexible bubble wrap, cut the toes out of tube socks and bunch them around the tubing.
3. Warm the air in the room.
4. Turn the temperature control on your HC220 Series down slightly. Gradually decrease the setting until the condensation no longer occurs.

Do I have to use distilled water with my HC220 Series?

The use of distilled water will maximize the life of the HC345 water chamber. Water from the faucet will often contain minerals which can damage the chamber, causing pitting in the base, even if the water is passed through a filter.


Can I use my HC220 Series in other countries?

Yes. Simply use the appropriate electrical socket adaptor and the HC220 Series will automatically adapt to any voltage supply.

Note: When moving your HC220 Series, ensure that the water chamber is empty. Machine failure due to water damage is not covered by warranty.

How do I replace the air filter?

Remove the filter from the holder as shown:

Note: When inserting the new filter, ensure that the black dot is facing towards the machine.

Can I use oxygen with my HC220 Series?

Yes, oxygen can be administered at the mask. Turn the HC220 Series on before turning on the oxygen. Ensure that the oxygen source is turned off before turning off the HC220 Series to avoid oxygen accumulating in the machine.

When do I replace my HC345 water chamber?

Replace the chamber if the plastic walls of the chamber become cracked or cloudy or the chamber base becomes pitted. (Note: If the chamber leaks at all, it should be replaced immediately).

| Product Specifications | | | |
|--|--|---|--|
| Dimensions: | 5.3" x 6.5" x 11.6" (135mm x 165mm x 295mm) | Temperature Control Settings of 1 – 4 Heater Plate Temperature operating maximum of 149°F (65°C) | |
| Weight: | 5 lbs (2.3kg) | | |
| Storage/transport temperature: | 14°F to 140°F (-10°C to 60°C) | | |
| Electrical Ratings | | Performance | |
| Supply Frequency: | 50-60 Hz | Delivered pressure: | 3-18cmH ₂ O |
| Supply Voltage and Current: | 100V-115V~, 1.05A 220V-240V~, 0.60A Note: These values represent average current | Max. delivered gas temperature: | 90°F (32°C) |
| | | Delivered Humidity: | 0-23mg/L |
| | | Warm-up time: | 20 minutes |
| | | Standards Compliance | |
| Suitable for use with a DC to AC inverter: | 200W minimum | Designed to conform to the requirements of: | IEC 60601-1 UL1431 CSA C22.2/No. 68 AS3200.1.0 EN60601-1 |
| Heater Plate: | 85W | | |

Troubleshooting

In the unlikely event that you tip water into the blower outlet:


- Remove air filter and chamber
- Turn unit on and run at prescribed pressure for 5 minutes
- Re-assemble and use as directed

If you feel that your HC220 Series is not operating correctly, consult your homecare dealer. Do not attempt to repair the HC220 Series yourself.

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